



www.MyBloodFirst.com

A METHOD TO REDUCE ALLOGENEIC BLOOD EXPOSURE AND HOSPITAL COSTS WHILE PRESERVING CLOTTING FACTOR CONCENTRATION AFTER CARDIOPULMONARY BYPASS

Keith A. Samolyk CCP, LCP, Global Blood Resources LLC, Somers, CT USA 06071 and Scott Beckmann CCP, Salem Hospital, Salem, OR. USA



Global Blood Resources LLC

Introduction

Cardiovascular Surgery remains responsible for as much as 20% of all transfusions in the United States despite recent data demonstrating that transfusions are independently linked to increased short and long term morbidity and mortality. (1, 12-13)

ECC circuits have long been viewed as a contributor to hemodilution. Condensed circuitry with prime volumes of 1,000-1,500 ml are now the norm and can be RAP'd (Retrograde Auto Primed) to reduce the hemodilution even further. (15) Blood volume remaining in the ECC at aortic decannulation has been traditionally salvaged by either processing with a "cell saver" or "chasing the ECC volume into the patient." (2-5, 7)

Cell processing conserves RBC's but discards plasma proteins. (8-10) Chasing the pump contents into transfer bags for infusion or directly into the patient stresses the kidneys to process extra fluid in a patient that is already volume overloaded. This stress may contribute to further organ dysfunction compared to maintaining normovolemic homeostasis. (11-12)

Observational data and descriptive statistics from a case series is presented to illustrate the use of the Hemobag® system. [See Figure One]

Method

A new blood conservation method and technology for blood salvaging, the Hemobag® (Global Blood Resources, Somers, CT 06071) deals directly with ECC volume at aortic decannulation.

The Hemobag® recovers and concentrates essentially all autologous whole blood, platelets and proteins from the ECC in a timely fashion for infusion, while maintaining the integrity and security of a safe primed circuit at all times. Use of the Hemobag® circuit allows for conventional ECC ultrafiltration during the procedure and works with any commercial hemofiltration.

After IRB approval a total of 66 patients undergoing cardiac surgery with CPB at Salem Hospital (Salem, Oregon) were randomly selected to the use of the Hemobag® Blood Salvage Device. Dependent variables and outcomes were compared between the Hemobag® treatment group and the control group.


Figure Two explains the method in more detail.

Table One

Parameter	Control Group	Hemobag® Group	p Value
Patient group size	66	66	NS
Percent male	70	74	NS
Age in years	66 +/- 11	66 +/- 13	0.300
BSA m ²	1.92 +/- 0.23	2.02 +/- 0.24	0.232
Pre-op weight kg	81 +/- 16	89 +/- 18	0.470
% Distribution of 5 surgeons	61/37/0/2/0	56/37/1/5/1	NS
% CABG surgery patients	66	65	NS
% Valve surgery patients	18	18	NS
% Valve + CABG surgery patients	16	17	NS
National Bayes risk score	5.4 +/- 6.7	5.0 +/- 6.6	0.782
LV ejection fraction %	59 +/- 16	56 +/- 22	0.456
CPB time min	131 +/- 60	131 +/- 42	0.989
Ischemic min	91 +/- 39	91 +/- 31	0.994
Pre-op HCT %	40 +/- 5	40 +/- 5	0.691
Pre-op platelet K/mm ³	227 +/- 81	224 +/- 89	0.809

Hemobag® Whole Blood "Recovery Loop"

Figure One



TS3 Tubing Set "Standard Loop"



The HEMOBAG® Blood Salvage Device is a reservoir system that allows the patient's whole blood to be Salvaged, Hemoconcentrated and Infused back to the same patient quickly, safely and efficiently in the same convenient reservoir bag (Insuring ECC integrity).

Hemobag® Case Series

Salem Hospital, Salem OR

- The end-CPB circuit blood for a group of 66 patients was processed using the Hemobag® (HB) device and technique
 - HB procedures and patients were selected randomly from all comers
- The Hemobag® patient results were compared to a concurrent matched control group (NHB)
 - The control group patients were selected to match the HB group patients by procedure, surgeon, CPB time, age and BSA
 - The control group end-CPB circuit blood was processed by the Cell Saver®
- Compared pre-op, operative, ICU and post-op parameters and outcomes between groups

Figure Two

Results

The average volume returned to the patient from the Hemobag® was 852 mL (1 SD = 197 mL).

The average time to hemoconcentrate the Hemobag® was just over 10 minutes.

Total blood protein concentration in the Hemobag® contents was 8.2 +/- 1.9 gm/dL.

Results are included in Tables One & Two and Figures Three - Five.

Hemobag® Case Series

Salem Hospital, Salem OR

HB vs. NHB group comparison:

- No significant difference in distribution of surgeons, procedures, patient age, BSA or gender between groups
- No significant difference in CPB and ischemic (clamp) times between groups
- No significant difference in National Bayes Risk Scores
- Same volume of pre-CPB autologous blood drawn (ANH)
- Same 100% use of cell-processor in both groups except the use of the Hemobag® to process the end-CPB circuit volume in the HB group

Figure Three

Hemobag® Case Series

Salem Hospital, Salem OR

HB vs. NHB group comparison:

- There were no significant difference in operative hematocrit nadir, post-op HCT or chest tube drainage volume
- Though not statistically significant there were strong tendencies toward clinically significant differences with sizeable financial consequences:
 - More Hemobag® patients received no blood products
 - Hemobag® patients received about 1/3 less total donor exposures on the average compared to control patients
 - Average RBC transfusion rates per patient were lower (50%) in the HB group (p = 0.053) compared to the control patients
 - Hemobag® patients blood costs were about 50% of control group blood costs per patient (p = 0.058): \$41,650. difference

Figure Four

Parameter	Control Group	Hemobag® Group	p Value
Hemobag® processed cc	NA	852 +/- 197	NA
Post-op bleeding cc/kg	9.4 +/- 6.5	7.9 +/- 7.2	0.191
FFP units per patient	1.1 +/- 2.4	0.9 +/- 1.5	0.546
Platelet pheresis units per patient	0.7 +/- 1.6	0.5 +/- 0.8	0.410
RBC transfusions per patient	1.8 +/- 2.9	1.1 +/- 1.6	0.053
Donor exposures per patient*	4.3 +/- 8.1	2.5 +/- 3.4	0.097
% Patients with no transfusions	41	50	
Cost blood products \$ per patient	1,417 +/- 2,487	777 +/- 1,055	0.008
Cost of blood products \$ for group	92,166	50,516	NA
Discharge % HCT	32 +/- 4	31 +/- 7	0.523
Pre-op creatinine mg/dl	1.4 +/- 2.2	1.1 +/- 0.4	0.222
Post-op creatinine mg/dl	1.6 +/- 1.0	1.8 +/- 1.9	0.524
Weight kg in ICU	NR	95.0 +/- 19.5	NR
ICU - pre-op weight change kg	NR	4.2 +/- 2.7	NR
Ventilated hours**	33 +/- 80	13 +/- 15	0.053
ICU hours**	72 +/- 98	53 +/- 53	0.161
Total hospital days**	11 +/- 12	9 +/- 4	0.163

Numerical data evaluated by chi-square analysis; other data analyzed by ANCOVA. *1 and NS are not significant at p < 0.05. NR is not measured, NR is not recorded and NA is not applicable. **One outlier (44 +/- 235) in each treatment group was removed before analysis. ***Three outliers were removed from each treatment group before analysis.

Table Two

Parameter	Control Group	Hemobag® Group	p Value
Pre-op HCT %	40 +/- 5	40 +/- 5	0.991
Pre-op platelet K/mm ³	227 +/- 81	224 +/- 89	0.809
Hemobag® content platelet K/mm ³	NR	238 +/- 73	NR
Post-CPB platelet K/mm ³	NR	121 +/- 46	NR
Post-op platelet K/mm ³	93 +/- 30	108 +/- 43	0.022
% of baseline post-op platelet count	55 +/- 15	51 +/- 16	0.079
Hemobag® content fibrinogen mg/dl	NR	461 +/- 174	NA
Post-CPB fibrinogen mg/dl	NR	206 +/- 69	NA
Low CPB %	31 +/- 5	31 +/- 7	0.796
Pre-CPB autologous blood draw cc/kg	4.9 +/- 2.4	4.3 +/- 2.3	0.594
Total heparin dose K IU/kg	NR	831 +/- 267	0.991
Hemobag® content HCT %	NA	44 +/- 6	NA
Low operative HCT %	23 +/- 3	24 +/- 3	0.031
% of baseline drop to low HCT	-41 +/- 10	-40 +/- 9	0.270

Numerical data evaluated by chi-square analysis; other data analyzed by ANCOVA. *1 and NS are not significant at p < 0.05. NR is not measured, NR is not recorded and NA is not applicable.

Discussion

Infusion of the CPB circuit residual blood concentrate appears to safely recover proteins, clotting factor and cell volume for all types of cardiac procedures which leads to reduced patient donor exposures, improved outcomes and reduction in the related costs.

Use of the Hemobag® allowed the clinicians to capture blood platelets and proteins that would have been normally discarded. (2, 5, 8-9) **Factor VII levels in three Hemobag® contents averaged a 259% increase.** (14)

The Hemobag® offers a new way to safely and efficiently manage and salvage autologous extracorporeal (ECC) whole blood for patients. Use of this new technique offers advantages over the current technologies of salvaging blood from ECCs while offering the potential to improve patient outcomes. (2, 5-10)

Prospective clinical studies are being conducted to assess the advantages in patient outcomes and potential reduction of allogeneic blood product use during cardiac surgery with the Hemobag® Blood Salvage Device end-CPB technique.

Total cost of blood products was reduced in the HemoBag® group

All blood products for the 33 of 66 HemoBag® patients receiving blood cost about \$50,516 versus the \$92,166 for the 39 of 66 control group patients who received a blood product **\$41,650. Difference!**




Figure Five

Reference

- Green J, Reynolds P, Spiess B, Levin J, Sutherland M, Aron T, McCarthy H, DeAnda A, Kasrijan V. Blood conservation is safe and effective for primary coronary artery bypass grafting. *Anesth Analg*. 2004;98:SCA11-134.
- Boldt J, Zickmann B, Czeke A, et al. Blood conservation techniques and platelet function in cardiac surgery. *Anesthesiology*. 1991; 75(3): 426-32.
- Kiziltepe U, Uysal A, Corapoglu T, et al. Effects of combined conventional and modified ultrafiltration in adult patients. *Ann Thorac Surg*. 2001; 71(12): 684-93.
- Boga M, Islemoglu, Badak I, et al. The effects of modified hemofiltration on inflammatory mediators and cardiac performance in coronary bypass surgery grafting. *Perfusion*. 2000; 15(2): 143-50.
- Leyh RG, Bartels C, Joubert-Hubner E, et al. Influence of modified ultrafiltration on coagulation, fibrinolysis and blood loss in adult cardiac surgery. *Euro J Cardiothorac Surg*. 2001; 19(2): 145-51.
- Luciani GB, Metton T, Vecchi B, et al. Modified ultrafiltration reduces morbidity after adult cardiac operations: a prospective, randomized clinical trial. *Circulation*. 2001;104(12 Suppl 1): I253-9.
- Nakamura Y, Masuda M, Tushima Y, et al. Comparative study of cell saver and ultrafiltration nontransfusion in cardiac surgery. *Ann Thorac Surg*. 1990; 49(6): 973-8.
- Tanemoto K, Hamanaka S, Morita I, Masaki H. Platelet activity of residual blood remaining in the Cardiopulmonary bypass circuit after cardiac surgery. *J Cardiovasc Surg (Torino)*. 2004; Feb;45(1):27-30.
- Guo XY, Duan H, Wang JJ, Luo AL, Ye TH, Huang YG, et al. Effect of intraoperative cell saver use on blood sparing and its impact on coagulation function. *Zhongguo Yi Xue Ke Xue Yuan Xue Bao*. 2004; Apr;26(2):186-91.
- Sodrakyan A, Gondak K, Patel D, Elefteriades JA. Volume expansion with albumin decreases mortality after coronary artery bypass graft surgery. *Chest* 2003; Jun 123(6):1853-7.
- Karkouti K, Beattie S, Wijesundhara D, Chan C, Rao V, Dattilo K, Djajani G, Ivanov J, Karski J. The degree of hemodilution during cardiopulmonary bypass is related to renal failure in adult cardiac surgery. *Anesth Analg* 2004; 98:SCA1-134.
- Habbib RH, Zacharias A, Schwann TA, Rordan CJ, Durham SJ, Shah A. Adverse effects of low hematocrit during cardiopulmonary bypass in the adult: should current practice be changed? *J Thorac Cardiovasc Surg*. 2003 Jun;125(6):1438-50.
- DeFoe GR, Ross CS, Olmstead EM, Surgenor SD, Fillingim MP, Groom RC, Forest RJ, Pieroni JW, Warren CS, Boggsian ME, Krumholz CF, Clark C, Clough RA, Weidner PW, Lacey SJ, Leavitt BJ, Mannin CA, Chaisworth DC, Marshall P, O'Connor GT. Lowest hematocrit on bypass and adverse outcomes associated with coronary artery bypass grafting. Northern New England Cardiovascular Disease Study Group. *Ann Thorac Surg*. 2001 Mar;71(3):769-76.
- Karkouti K, Beattie WS, Wijesundhara DN, Yau TM, McCluskey SA, Ghanam M, Sutton D, van Rensburg A, Karski J. Recombinant factor VIIa for intractable blood loss after cardiac surgery: a propensity score-matched case-control analysis. *Transfusion*. 2005 Jun;45(1):26-34.
- Zelinka ES, Ryan P, McDonald J, Larson J. Retrograde autologous prime with shortened bypass circuits decreases blood transfusion in high-risk coronary artery surgery patients. *J Extra-Corpor Technol*. 2004 Dec;36(4):343-7.